DUROSIL L

Putty C-SILICONE IMPRESSION MATERIAL
Putty C-SILIKON ABFORMMATERIAL

DUROSIL S

Light Body C-SILICONE IMPRESSION MATERIAL
Dünnfließendes C-SILIKON ABFORMMATERIAL

PASTE HARDENER

CATALYST FOR C-SILICONE IMPRESSION
KATALYSATOR FÜR C-SILIKON ABFORMMATERIAL

EN

DE

IT


PRESIDENT DENTAL
Zehentsiedelweg 7 81247 München / Germany
Tel: 089-127 660 240 Fax: 089-127 660 269
www.presidentdental.com info@presidentdental.com

EN

Durosil L: Polysiloxane for the first impression characterised by initial flow in the working phase and high rigidity after hardening. Recommended for the stage technique (double impression)

Applications for soft and hard:

- Initial Impression for the double impression technique
- Situation impression - crush bite
- Spoon material for the sandwich technique


Dosage: Use the measuring spoon to dose the required quantity of Durosil L and spread it out over your hand (note: for the correct dosage the putty should be level with the top the measuring scoop). Press the rim of the measuring scoop onto the putty once for each measure dosed. For each measure of material, apply two strips of President Paste Hardener as long the measuring scoop (4 cm.)

Impression preservation: Store at room temperature under dry conditions

Storage: Store at not above 27°C

Mixing Time	: 0:30 min
Total Working Time	: 1:15 min
Time in Mouth	: 3:15 min
Setting Time	: 4:30 min
Strain Compression	: 2-5%
Recovery From Deformation	: >98%
Linear Dimensional Change After 24 hours	: <0,2%

Durosil S: Low viscosity, extremely fluid, hydrocompatible, polysiloxane, Recommended for the two stage technique (double impression)

Preparation and processing

Storage: Store at not above 27°C

Mixing Time	: 0:30 min
Total Working Time	: 1:15 min
Time in Mouth	: 3:30 min
Setting Time	: 5:00 min
Compressive deformation (min-max)	: 2-8%
Elastic recovery	: >99%
Dimensional stability	: <0,7%

Paste Hardener: Paste hardener for the condensation curing Durosil products Paste Hardener comes for silicone impression materials Able to be metered sparingly Can be applied well Easy to mix in

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